## Sybr Safe DNA Gel Stain

### Description

SYBR Safe DNA gel stain has been specifically designed for reduced mutagenicity, making it safer to replace the highly toxic ethidium bromide (EtBr) for staining DNA in agarose or acrylamide gels. SYBR Safe DNA gel stain showed no or very low mutagenic activity at concentrations well above working concentrations used for gel staining. Sybr Safe Nucleic Acid Gel Stain is highly sensitive than EtBr either as precast gel stains or post gel stains.

DNA bands stained with SYBR Safe DNA gel stain can be detected using a standard UV transilluminator, a visible-light transilluminator, or a laser-based scanner. The stain is also suitable for staining RNA in gels.

Sybr Safe Nucleic Acid Gel Stain 10,000X in DMSO is a concentrated Sybr Safe solution that can be diluted 10,000 times for use in precast gel staining for ~3,300 times for use in post gel staining according to the procedures described below. One vial(0.5ml) of 10,000X solution can be used to prepare at 100 precast minigels or post-stain at least 100 minigels.

Gel staining with Sybr Safe is compatible with downstream applications such as gel extraction and cloning. Sybr Safe is efficiently removed from DNA by phenol/chloroform extraction and ethanol precipitation.

Bound to nucleic acids, Sybr Safe DNA Gel stain has fluorescence excitation maxima at 280 and 502 nm, and an emission maximum at 530 nm.

### Storage

Store the SYBR Safe DNA gel stain at any temperature between 2 to 25°C. It must be completely thawed and mixed if the SYBR Safe is frozen at low temperatures before using. Repeated freeze-thawing has minimal impact on product performance. SYBR Safe DNA gel stain showed no or very low mutagenic activity.

### Features

- Safety: Nearly nonmutagenic and noncytotoxic
- Easy disposal: Safe to dispose in the drain
- Compatibility: Spectrally compatible with existing instruments
- Sensitivity: Higher signal but lower background
- Stability: can be stored at RT and microwavable

### Post-staining Protocol

- 1) Run gels as usual according to your standard protocol.
- 2) Dilute the Sybr Safe 10,000X stock reagent 3,300 fold to make a 3X staining solution TAE or TBE buffer prior to use.
- 3) Carefully place the gel in a suitable polypropylene container. Gently add a sufficient amount of the 3X staining solution to submerge the gel.
- 4) Agitate the gel gently at room temperature for ~30 minutes.
- 5) Image the stained gel with a 254nm transilluminator, a Dark Reader or a similar transilluminator, or a laser-based gel scanner using a long path green filter such as a SYBR filter or GelStar filter.
- 6) Staining solution can be reused at least 2~3times. Store staining solution at room temperature protected from light.
- 7) Stained gels can be viewed using a standard 300 nm transilluminator, a 254 nm epi- or transilluminator, or a blue-light transilluminator.

### Pre-cast Protocol

- 1) Prepare molten agarose gel solution using your standard protocol.
- 2) Dilute the Sybr Safe 10,000X stock reagent into the molten agarose gel solution at 1:10,000 and mix thoroughly. Sybr Safe can be added while the gel solution is still warm (50~60°C). For example if you run TBE gels and require 30 mL of molten agarose for your tray, add 3 µL of 10,000X SYBR Safe stain concentrate to 30 mL of 1X TBE.
- 3) Cast the gel and allow it to solidify. Any leftover gel solution may be stored and reheated later for additional gel casting. Sybr Safe precast gels may be stored at 4°C for later use.
- 4) Load samples and run the gels using your standard protocol.
- 5) Image the stained gel with a 254nm transilluminator, a Dark Reader or a similar transilluminator or a laser-based gel scanner using a long path green filter such as a SYBR filter or GelStar filter.
- 6) Stained gels can be viewed using a standard 300 nm transilluminator, a 254 nm epi- or transilluminator, or a blue-light transilluminator.

Note: The pre-cast protocol is not recommended for polyacrylamide gels.

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