
## Fluorescent Dye Isothiocyanates

Isothiocyanates form thioureas upon reaction with amines. It is proven that some thiourea products (in particular, the conjugates from  $\alpha$ -amino acids/peptides/proteins) are much less stable than the conjugates that are prepared from the corresponding succinimidyl esters. It has been reported that antibody conjugates prepared from fluorescein isothiocyanates deteriorate over time. We strongly recommend that you use succinimidyl esters for your conjugations whenever possible. There are few factors that need be considered when SE compounds are used for conjugation reaction: 1). Solvents: For the most part, reactive dyes are hydrophobic molecules and should be dissolved in anhydrous dimethylformamide (DMF) or dimethylsulfoxide (DMSO). 2). Reaction pH: The labeling reactions of amines with isothiocyanates are strongly pH dependent. Isothiocyanate reagents react with non-protonated aliphatic amine groups, including the terminal amines of proteins and the  $\epsilon$ -amino groups of lysines. Protein modifications by isothiocyanates may require a pH 9.0-10.0 for optimal conjugations. 3). Reaction Buffers: Buffers that contain free amines such as Tris and glycine must be avoided when using an amine-reactive reagent. Ammonium salts (such as ammonium sulfate and ammonium acetate) that are widely used for protein precipitation must also be removed before performing dye conjugations. High concentrations of nucleophilic thiol compounds should also be avoided because they may react with the labeling reagent to form unstable intermediates that could destroy the reactive dye. 4). Reaction Temperature: Isothiocyanate conjugations are usually done at room temperature. However, either elevated or reduced temperature may be required for a particular labeling reaction.

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
#### Features and Biological Applications


5(6)-TRITC is an amino-reactive labeling reagent that is widely used in preparing bioconjugates of proteins and nucleic acids. The resultant conjugates have similar spectral properties to those prepared from 5(6)-TAMRA, SE. However, the latter conjugates are much more stable. *Cautions must be exercised for the storage of FITC conjugates.*

#### References

1. Pellestor, F., *et al.*, Fast multicolor primed in situ protocol for chromosome identification in isolated cells may be used for human oocytes and polar bodies. *Fertil Steril* 2004, **81**, 408-15.
2. Gustafsson, M.K., *et al.*, No nerves and their targets in a tapeworm: An immunocytochemical study of cgmp in *hymenolepis diminuta*. *Parasitol Res* 2003, **90**, 148-52.
3. Takeno, S., *et al.*, Increased nitric oxide production in nasal epithelial cells from allergic patients--rt-pcr analysis and direct imaging by a fluorescence indicator: Daf-2 da. *Clin Exp Allergy* 2001, **31**, 881-8.
4. Meadows, D.L., *et al.*, Determining the extent of labeling for tetramethylrhodamine protein conjugates. *J Immunol Methods* 1991, **143**, 263-72.

### 5-TRITC [Tetramethylrhodamine-5 isothiocyanate]

Cat#	Size	Price	MW	Abs	Em	Soluble in	Storage
415	5 mg	\$145	443.52	543 nm	571 nm	DMF or DMSO	-20 °C and desiccated

#### Features and Biological Applications

5-TRITC (also called G isomer) is a purified single isomer of the 5(6)-TRITC mixed isomers. This labeling reagent is predominantly used in labeling peptides and proteins. *Cautions must be exercised for the storage of TRITC conjugates.*

#### References

1. Pellestor, F., *et al.*, Fast multicolor primed in situ protocol for chromosome identification in isolated cells may be used for human oocytes and polar bodies. *Fertil Steril* 2004, **81**, 408-15.
2. Takeno, S., *et al.*, Increased nitric oxide production in nasal epithelial cells from allergic patients--rt-pcr analysis and direct imaging by a fluorescence indicator: Daf-2 da. *Clin Exp Allergy* 2001, **31**, 881-8.
3. Newkirk RF and Mack J (1992). Improved indirect fluorescence immunocytochemical method using counter stains. *Biotechniques* **13**, 536-8.

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
## 6-TRITC; R isomer [Tetramethylrhodamine-6-isothiocyanate]

Cat#	Size	Price	MW	Abs	Em	Soluble in	Storage
417	5 mg	\$145	443.52	544 nm	572 nm	DMF or DMSO	-20°C and desiccated

### Features and Biological Applications

6-TRITC (also called R isomer) is the other isomer of the TRITC labeling reagent that is widely used in preparing bioconjugates of proteins and nucleic acids. Complimentary to 5-TRITC, the 6-isomer is predominantly used in labeling nucleotides and nucleic acids.

*Cautions must be exercised for the storage of TRITC conjugates.*

### References

- Kahn E, *et al.* (1999). Confocal-multilabeling, ultrasensitive TUNEL analysis of DNA breaks in individual cells. *Anal Quant Cytol Histol* **21**, 1-7.
- Nederlof PM, *et al.* (1992). Fluorescence ratio measurements of double-labeled probes for multiple in situ hybridization by digital imaging microscopy. *Cytometry* **13**, 839-45.
- Meadows, D.L., *et al.*, Determining the extent of labeling for tetramethylrhodamine protein conjugates. *J Immunol Methods* 1991, **143**, 263-72.

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Beijing Fluorescence Biotechnology Co. Ltd