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SafeRed[®] Nucleic Acid Gel Stain

SafeRed[®], 10000X GelRed[®], SafeRed[®]和DuGreen[®] are registered trademarks and strictly forbidden to be used without authorization.

Product Description

SafeRed[®] Nucleic Acid Gel Stain is the third generation nonmutagenic fluorescent nucleic acid gel stain to replace the highly toxic traditional ethidium bromide (EtBr) for staining dsDNA, ssDNA or RNA in agarose gels or polyacrylamide gels. The third party Ames test by WuXi Pharmacy shows that SafeRed and DuGreen are nonmutagenic for gel staining. DuGreen and SafeRed Nucleic Acid Gel Stain are highly sensitive than EtBr either as precast gel stains or post gel stains.

We can directly replace EB with SafeRed[®] without changing the existing imaging system since SafeRed[®] and EB have the similar spectra. SafeRed[®] can also be used to stain dsDNA, ssDNA or RNA in polyacrylamide gel via post gel staining. Gel staining with SafeRed[®] is compatible with downstream applications such as gel extraction and cloning. SafeRed[®] is efficiently removed from DNA by phenol/chloroform extraction and ethanol precipitation.



Features

- Safety: Nonmutagenic and noncytotoxic
- Easy disposal: Safe to dispose in the drain
- Compatibility: Spectrally compatible with all the existing instruments
- Sensitivity: Higher signal but lower background
- Stability: can be stored at RT and microwavable

Materials needed:

Agarose, TAE Buffer, 6X Sample Loading Buffer, DNA ladder standard, Electrophoresis chamber, Power supply, Gel casting tray and combs, Staining tray, Gloves, Pipette and tips

Precast Protocol for Agarose Gel Electrophoresis

Pouring a Standard 1% Agarose Gel:

- 1 Pour 1g of Agarose powder and 100mL of 1xTAE or 1xTBE into a glass flask.
Note: Agarose gels are commonly used in concentrations of 0.7% to 2% depending on the size of DNA bands.
- 2 Melt the agarose in a microwave for 1-3min (until the agarose is completely dissolved and the solution becomes clear, the solution is completely clear and no small floating particles are visible.)
Note: Caution HOT! Be careful! Do not let the solution boil for long periods as it may boil out of the flask.
- 3 Let agarose solution cool down to about 50-55°C, Add 10uL SafeRed per 100mL gel, gently stirring the gel solution to mix thoroughly.
- 4 Set horizontal gel apparatus. Pour molten agarose onto gel plate to a depth of 4~8 mm. Insert a comb until its base is 1 mm from the base of the gel. Allow cooling.
Note: Pour slowly to avoid bubbles, which will disrupt the gel. Any bubbles can be pushed away from the well comb or towards the sides/edges of the gel with a pipette tip.
- 5 Place newly poured gel at 4°C for 10-15 minutes OR let sit at room temperature for 20-30 minutes, until it has completely solidified. Remove comb and submerge in 4X TAE buffer

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Loading Samples and Running an Agarose Gel:

- 1 Once solidified, place the agarose gel into the gel box.
- 2 Add 1/5 volume 6X Loading buffer to sample.
Note: Loading buffer serves two purposes: 1) it provides a visible dye that helps with gel loading and will also allow you to gauge how far the gel has run while you are running your gel; and 2) it contains a high % glycerol, so after adding it your sample is heavier than water and will settle to the bottom of the gel well, instead of diffusing in the buffer.
- 3 Fill gel box with 1xTAE (or TBE) until the gel is covered.
- 4 Carefully Pipette 5 μ l of the DNA ladder standard into the first lane of the gel.
Note: When loading the sample in the well, maintain positive pressure on the sample to prevent bubbles or buffer from entering the tip. Place the very top of the tip of the pipette into the buffer just above the well. Very slowly and steadily push the sample out and watch as the sample fills the well. After all samples are unloaded, push the pipettor to the second stop and carefully raising the pipette straight out of the buffer.
- 5 Carefully pipette 5 μ l of each sample/Sample Loading Buffer mixture into other separate wells in the gel.
- 6 Run the gel at 80-150V until Bromophenol Blue is near the end of the gel. This dye runs at ~800b.
Note: Black is negative, red is positive. (The DNA is negatively charged and will run towards the positive electrode.) **Note:** A typical run time is about 0.5-1 hour, depending on the gel concentration and voltage.
- 7 Turn OFF power, Disconnect the electrodes from the power source, and then carefully remove the gel from the gel box.
- 8 Visualize or Image the stained gel with a standard transilluminator (302 or 312nm), and photograph the gel using an ethidium bromide filter.

Post-Staining Protocol

1. Run Agarose Gel Electrophoresis according to hereinbefore.
2. Dilute 10000X SafeRed stock solution to 3X staining solution in H₂O with 0.1M NaCl(e.g, add 15ul of SafeRed[®] 10,000X stock reagent and 5ml 1M NaCl to 45ml H₂O). Note: including 0.1M NaCl in the staining solution enhances sensitivity, but may promote dye precipitation if the gel stain is reused. (This solution can be reused at least 2~3 times protected from light).
3. Using gloves, carefully remove the gel from the casting tray and place into the staining dish, add 3X SafeRed staining solution to submerge the gel.
4. Agitate the gel gently at room temperature for~30 minutes and the amount of the stain may depend on the thickness of the gel and the percentage of the agarose with gently shaking.
5. Rinsing the gel with water can reduce the background.
6. Visualize or Image the stained gel with a standard transilluminator (302 or 312nm), and photograph the gel using an ethidium bromide filter.

Tips and FAQ

◆ How do you get better resolution of bands?

A few simple ways to increase the resolution (crispness) of your DNA bands include: a) running the gel at a lower voltage for a longer period of time; b) using a wider gel comb; or c) loading less DNA into well.

◆ How do you get better separation of bands?

If you have similarly sized bands that are running too close together you can adjust the agarose percentage of the gel to get better separation. A higher percentage agarose gel will help resolve smaller bands from each other, and a lower percentage gel will help separate larger bands.

◆ 10% Rule:

For each sample you want to load on a gel, make 10% more volume than needed because several microliters can be lost in pipetting. For example, if you want to load 1.0 μ g in 10 μ L, make 1.1 μ g in 11 μ L.

◆ Diffused strip or smeared DNA band(s).

SafeRed and DuGreen cannot go into the DNA double helix structure like EB because their big molecular weight is more than 1000. Please try the following methods to reduce the smeared DNA band(s) or diffused strip: 1) Load half of your DNA sample and DNA marker; 2) Pour a lower percentage agarose gel for better resolution of large fragments; 3) Using the fresh TAE buffer or TBE buffer instead of TAE. 4) Please run the gel stain by the post-staining protocol.

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Beijing Fluorescence Biotechnology Co. Ltd