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DRAQ5[®] Fluorescent Probe

DRAQ5[®]和 DRAQ7[®] are registered trademarks and strictly forbidden to be used without authorization.

DRAQ5[®] Fluorescent Probe, 5mM, 50 μ L, sufficient for ~4 x 96 wells

DRAQ5[®] Fluorescent Probe, 5mM, 200 μ L, sufficient for ~16 x 96 wells

MW~413 Da, Wavelength: Ex/Em= 646/697nm

Features

- Stable at room temperature
- Rapid staining of cells without wash steps
- No UV excitation and no emission spectral overlap with PE

Storage&Shipping:

Store at 4°C protected from light. Product is shipped at ambient temperature.

Discription:

The DRAQ5[®] Fluorescent Probe is a far-red DNA stain for use in live or fixed cells. Because of its far-red excitation and emission, the DRAQ5 Fluorescent Probe can be multiplexed with many other fluorophores and is ideal for cells expressing green fluorescent protein (GFP) fusion proteins. DRAQ5 Fluorescent Probe is compatible with many existing protocols across a wide range of instrumentation platforms. This fluorescent probe can be used for cellular imaging.

Fluorescence cell-based assays, such as flow cytometry, in-cell ELISA, fluorescence microscopy and high-content imaging require a fluorescent label to identify individual cells. When using multiple fluorescent probes to detect different cellular targets or activities, each probe must have a fluorescent spectrum different than the other probes. The blue-fluorescent DNA-binding probes, Hoechst and DAPI are frequently used; however, these probes cannot be used when UV illumination is unavailable or other blue-emitting fluorescent probes are used. Therefore, nuclear probes that emit in a color other than blue are useful for cell identification and counting, and for determining nuclear morphology and DNA content.

The DRAQ5[®] Fluorescent Probe emits in the far-red region, is lipophilic and crosses cell and nuclear membranes in live and fixed cells and tissues for rapid DNA staining. This stain is water-soluble, supplied ready to use and does not require RNase, cell lysis, or a washing step, making it compatible with automation. Because DNA staining is stoichiometric, the DRAQ5 Fluorescent Probe can be used for DNA content analysis in cell proliferation studies.

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Procedure for DNA Staining

Note: DRAQ5® Fluorescent Probe is usually added as the last stain in a labeling procedure because no washing is required. Alternatively, add this probe in assay medium for a live cell assay. Read the supplied material safety data sheet before handling DRAQ5® Probe.

1. Prepare phosphate-buffered saline (PBS, without sodium azide) or the appropriate culture media for the specific cells.
2. Resuspend cells in PBS or media at $\leq 4 \times 10^5$ cells/mL in a test tube. For adherent cells estimate the number of cells based on confluence level or tissue section dimensions.
3. Add DRAQ5 Fluorescent Probe directly as supplied according to the volumes indicated in Table 1. Add directly on top of tissue sections and adherent cells or add DRAQ5 Fluorescent Probe in fresh media.
4. Gently mix and incubate for 5–30 minutes at room temperature. DRAQ5 Fluorescent Probe staining is accelerated at 37°C and may be reduced to 1–3 minutes. For time-lapsed assays (e.g., studying translocation of an EGFP-tagged protein) DRAQ5 Fluorescent Probe may be added to the assay media for the duration of the assay (typically 0.5–3 hours) at 1µM before adding any agonist/antagonist.

Note: Protect cells from light during incubation if other (immuno-) fluorescent stains have been applied to the cells before the DRAQ5 Fluorescent Probe staining.

5. Cells can be analyzed directly without further treatment or washing.

Cell sample preparation:		VOLUME OF DRAQ5® (AS SUPPLIED) REQUIRED FOR A CONCENTRATION OF:		
No. of cells:	in volume:	5 µM	10 µM	20 µM
1 x 10 ⁶	2500 µl	2.5 µl	5 µl	10 µl
4 x 10 ⁵	1000 µl	1 µl	2 µl	4 µl
2 x 10 ⁵	500 µl	0.5 µl	1 µl	2 µl
1 x 10 ⁵	250 µl	0.25 µl	0.5 µl	1 µl
5 x 10 ⁴	125 µl	0.13 µl	0.25 µl	0.5 µl

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Beijing Fluorescence Biotechnology Co. Ltd