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## DRAQ7<sup>®</sup> Fluorescent Probe

DRAQ5<sup>®</sup>和 DRAQ7<sup>®</sup> are registered trademarks and strictly forbidden to be used without authorization.

DRAQ7 Fluorescent Probe, 0.3mM, 1mL

MW~413 Da, Wavelength: Ex/Em= 646/697nm

### Features

- Stable at room temperature
- Rapid staining of dead cells without wash steps
- No UV excitation and no emission spectral overlap with PE

### Storage&Shipping:

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze. Product is shipped at ambient temperature.

### Discription:

DRAQ7<sup>®</sup> is a bright, easy-to-use far-red fluorescent DNA dye. DRAQ7<sup>®</sup> Dye is a membrane-impermeable dye that rapidly stains the double-stranded DNA (dsDNA) of dead or fixed cells for cell cycle analysis by DNA content, nuclear visualization, or discrimination of nucleated cells from debris or enucleated cells. Because DRAQ7<sup>®</sup> is impermeable to intact cells, it may also be used as a viability dye. DRAQ7<sup>®</sup> has excitation maxima at 599/644 nm and emission maxima at 678/697 nm, but can also be sub-optimally excited by the 488 nm. Its emission wavelength maximum is 678 nm, or 694 nm when intercalated with double-stranded DNA.

### Recommended Procedure for DNA Staining

#### 1. Staining of Live Cells for Viability Analysis by Flow Cytometry

1.1 Obtain a single cell suspension.

1.2 Resuspend cells in 1X Dulbecco's Phosphate Buffered Saline (DPBS) or other azide-free buffer containing 1-3 μM DRAQ7<sup>®</sup>.

- a. The optimal concentration of DRAQ7<sup>®</sup> for viability analysis may vary by cell type. We recommend titrating the reagent for your cell type of interest in early experiments.
- b. Additionally, apoptotic cells may stain with variable amounts of DRAQ7<sup>®</sup>. We recommend co-staining with, eg, our FITC Annexin V if further analysis of apoptotic cells is desired.

1.3 Incubate 5 minutes at room temperature. No wash is necessary prior to analysis.

1.4 Proceed to analysis by flow cytometry.

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## 2. Staining of Fixed Cells for DNA Content Analysis by Flow Cytometry

- 2.1. Obtain a single cell suspension.
- 2.2. Treat cells on ice for 30 minutes with 70–80% ice–cold ethanol. Ethanol fixation typically provides the most resolved histograms.
- 2.3. Wash cells once with PBS.
- 2.4. Dilute DRAQ7<sup>®</sup> to 20  $\mu$ M in 1X DPBS or other azide–free buffer immediately prior to use.
- 2.5. Stain cells for 5–15 minutes at a cell density of 0.5E6 cells/mL or less. No further wash is necessary prior to analysis. The optimal cell density and concentration of DRAQ5<sup>®</sup> for DNA content analysis may vary by cell type. Assay conditions should be optimized in early experiments for best results.
6. Proceed to analysis by flow cytometry.

## 3. Immunofluorescent Staining of Fixed Cells for Nuclear Visualization

- 3.1. Fix and permeabilize cells as desired.
- 3.2. Dilute DRAQ7<sup>®</sup> solution to 5–20  $\mu$ M in 1X DPBS or other azide–free buffer immediately prior to use.
- 3.3. Add DRAQ7<sup>®</sup> solution to each sample at least 5 minutes before analysis.
- 3.4. Proceed to imaging. We recommend using a 715LP or longer wavelength filter, though the dye is well–detected in filters typically used to detect Super Fluor. 647 (eg, 660/20 or 692/40).  
Note that dsDNA–bound dye will fluoresce brightly in the nucleus and unbound dye may fluoresce dimly in the cytoplasm, allowing segmentation of the cytoplasmic and nuclear compartments.

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Beijing Fluorescence Biotechnology Co. Ltd